

## Upholstered Drop Arm Wheeled Commode

HCPCS: E0165



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Simple, tool-free assembly
- Fully knock-down for easy storage and transport
- Drop-arms allow for easier lateral transfers
- Swing-away detachable footrests
- Collection pan included secures underneath the seat and can be easily removed from the backside of the unit
- Comfortable, wide backrest
- Durable and easy-to-clean padded vinyl seat cover

PRODUCTS	Description
11120SV-1F	Wheeled Commode with Footrests, 1/CV
SPECIFICATIONS	Description
Carton Shipping Dimensions	22.25" (L) x 22.45" (W) x 24.02" (H)
Carton Shipping Weight	30.87 lbs
HCPCS	E0165
Product Weight	27.99 lbs
Seat Dimensions	16" (W) x 16" (D)
Warranty	Limited Lifetime
Weight Capacity	350 lbs
Width Between Arms	18.5"   Outside Arms: 22.25"
Material Composition	Steel