

Med-Aire 8" Alternating Pressure and Low Air Loss Mattress System

HCPCS: E0277



FIGURE A



FIGURE B

- Active mattress replacement system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Control unit offers 10 minute cycles when alternating and can also be used in static mode
- Adjustable patient weight settings allow for optimal immersion and patient comfort & compliance
- Audio and Visual alarm alerts when air pressure is low
- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
- 8" deep air cells are easily removed and replaced with air leakage plugs
- Removable Nylon stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS	Description
14027	Med Aire Low Air Loss, 8in 1/cs

SPECIFICATIONS	Description
Control Unit Airflow LPM	8 LPM
Control Unit Alarms	Visual / Audible
Control Unit Cycle Time	10 Min.
Control Unit Type	Compressor
Dimensions	36" (W) X 80" (L) X 8" (H)
Mattress Cover Material	Fluid-Resistant Nylon
Waranty	1 Year Limited
Weight Capacity	350 lbs

ACCESSORIES	Description
14027M	Mattress Only for 14027, 1/ea, EA
14427QF	Quick Fill Pump w/adapters, 1 EA/12 CV, EA

PARTS	Description
14027-BASE	Replacement Base for 14027, 1/ea, EA
14012	CPRPlugForLowAirLossMatt,2Port,1ea, EA
14007	Tubing Adapter/APP,1/ea, EA
14027P	Pump for 14027, 1/ea, EA
14001T	Tubing for 14001/14002,1/ea, EA

PARTS	Description
14027-CELL	Repl Bladder - 14027, 1/ea, EA
14027c	Cover for 14027,1/ea, EA
14027-cell1	Repl Bladder w/Lsr Holes - 14027, 1/ea, EA
14028v	OneWayValveForLowAirLossMatt,1/ea, EA
14028t	Tubing,7", LAL, 1/ea, EA
14027-fftc	Connector for 14027 Cell, Female, 1/ea, EA
14027-fmtc	Connector for 14027 Cell, Male, 1/ea, EA
14027-fbs	Stopper for 14027, 1/ea, EA
14028h	Hose only for 14028, 14027, 1/ea, EA
14028s	Pump Socket,2 Port, 14027,14028,1/ea, EA
14027-fk	Control Knob,14027 W/Ser #4E1506 & Newe, EA
140pdk	Pressure Dial Knob, LAL Pump, 1/ea, EA
14027-fph	Pump Hooks, 14027, 1/set, EA
14029f	Filter only for 14029,1/ea, EA
14530pf	Pump Filter, 14530P, 1/ea, EA