


AL560 / AL560XL

AUTOMATIC UNIVERSAL POWER WHEELCHAIR LIFT



ADJUSTABLE POWER CHAIR SOLUTION

The easy way to secure a wide range of popular power chairs. Hold-Down arm is padded to hold securely while cushioning chair. The 26.25" x 38.5" platform (AL560) and 30.75" x 38.5" (AL560XL) platform adjusts for a custom fit. Manual crank backup and license plate mount are included.

Swing-away option available.

AL560 PRODUCT SPECIFICATIONS

| | |
|----------------------------|--|
| LIFTING CAPACITY | 350 lb |
| ACCOMMODATES | Power wheelchairs |
| INSTALLED WEIGHT | 89 lb |
| PLATFORM SIZE | 26.25" x 38.5" adjustable platform |
| HITCH CLASS | II or greater |
| POWERED | Fully powered |
| WARRANTY | 3-year |
| INCLUDED | <ul style="list-style-type: none">- AL107 hitch post- AL123 class III hitch adapter- License plate mount |
| SECUREMENT | <ul style="list-style-type: none">- Padded hold-down arm- Wheel cradles |
| AVAILABLE OPTION(S) | <ul style="list-style-type: none">- Battery pack- Cover- Swing-away |

AL560XL PRODUCT SPECIFICATIONS

| | |
|----------------------------|---|
| LIFTING CAPACITY | 350 lb |
| ACCOMMODATES | Larger power wheelchairs |
| INSTALLED WEIGHT | 101 lb |
| PLATFORM SIZE | 30.75" x 38.5" adjustable platform |
| HITCH CLASS | III or greater |
| POWERED | Fully powered |
| WARRANTY | 3-year |
| INCLUDED | <ul style="list-style-type: none">- AL107 hitch post- AL123HD class III heavy duty hitch adapter- License plate mount |
| SECUREMENT | <ul style="list-style-type: none">- Padded hold-down arm- Wheel cradles |
| AVAILABLE OPTION(S) | <ul style="list-style-type: none">- Battery pack- Cover- Swing-away |



Harmar.com
800.833.0478

17SEP2020 | MKT-000176 REV A

