

**ALL-NEW**  
**Sporty Fold-Flat Scooter!**  
**Airline Friendly\*!**  
**USB Port Standard!**



Seat Color  
 Black Vinyl

Body Colors  
 Blue  
 Red

LIGHTWEIGHT PORTABLE SCOOTER



GB120

The Buzzaround CarryOn Fold-Flat Scooter is small in stature and big on fun! This super portable, sporty new scooter will take you where you want to go in comfort and style! It features an infinite adjustable tiller, full front and rear light package, lithium ion battery pack (travel battery pack sold separately), and 300 lb. weight capacity. The heaviest piece is 47.5 lbs. (without seat and batteries) and the height when folded is only 12". It disassembles easily into two pieces and is designed to store in the smallest compartments. Pack your dreams and go anywhere with the CarryOn! \*Optional 6.5 AH battery pack can be removed and taken as a carry on while traveling by air, sea or train.





LIGHTWEIGHT PORTABLE SCOOTER

<b>Model</b>	<b>GB120</b>
Medicare Code	Not Coded
Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	3.6 mph
Operating Range <sup>1</sup>	18 miles (15 AH) 7.8 miles (6.5 AH)
Ground Clearance Under Deck	4"
Ground Clearance to Transaxle/Motor	1.5"
Turning Radius	48.75"
Type Batteries (Lithium Ion)	15 AH (Standard) 6.5 AH (Optional)
Color	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	50 Amp
<b>Dimensions</b>	
Length	38"
Width	19.5"
Height (ground to top of back of seat)	34.5"
Maximum Height (folded without seat)	12"
Ground to Top of Seat	23.5"
Ground to Top of Deck	6"
Deck to Top of Seat	17.5"
Front Axle to Rear Axle	28.7"
<b>Weight of Unit: (assembled)</b>	65.5 lbs
Main Body	47.5 lbs
Seat w/Arms	12 lbs
Batteries	15 AH 6 lbs 6.5 AH 3.8 lbs
<b>Tires:</b>	Solid/Low Profile
Front	6"
Rear	7"
Rear Anti-Tip Wheels	2"
<b>Standard Seat</b>	Stadium Style
Back Height - no headrest	13.25"
Width x Depth	15" x 13.75"
Color	Black Vinyl
<b>Warranty:</b>	
Frame	Lifetime
Drive Train	2 Years
Electronics <sup>2</sup>	2 Years
Batteries	1 Year



Tiller & self-locking seat frame fold for a compact profile



Sporty new look for on the go!



Front & rear lighting package



Front carry handle makes travel easy



Conveniently located USB port



LED battery indicator on the console



Red & blue color shrouds available



Features

- Designed for minimal storage space
- 4-point carry handles for easy transport
- 12" height when folded (without seat)
- 4" ground clearance
- 300 lb. weight capacity
- Flip up armrests
- Infinite adjustable tiller
- USB port & water bottle holder on tiller
- Disassembles easily
- Long range coverage up to 18 miles
- Easy access 15 AH lithium ion battery pack standard
- Optional airline friendly 6.5 AH lithium ion battery pack removable for carry on (Available Late Summer 2020)

Order the optional Golden Travel Case for on the go storage!



Scan for Video

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

