

RevVA3

VISION SPORT

FULL SIZE



merits

8/14/2020

FEATURES

- True mid-wheel-drive: 20.6" turning radius
- Robust independent suspension can absorb up to 250 lbs. on the footplate without tipping forward
- Premium, semi-reclining ergonomic seat; up to 135°
- 10" knobby drive tires for superior traction
- Dynamic LiNX controller for convenient Bluetooth programming and diagnostics
- 16" seat width with 16"-18" depth adjustability
- 18" or 20" seat widths with 18"-20" depth adjustability
- Red and Blue color panels included

SPECIFICATIONS

Weight Capacity	300 lbs
Maximum Speed	5 mph
Per-Range Change	18 miles
Batteries	U1 - 34Ah
Ground Clearance	2"
Turning Radius	20.6"
Overall Length	34"
Overall Width	24"
Drive Wheels	10"
Front Casters	6"
Rear Casters	6"
Suspension	TruStep
Seat Width	16", 18" or 20"
Seat Depth	16", 17", 18", 19" or 20"
Backrest Height	18" or 20"
Power Seat Lift	P326D : 7.5"
Controller	P326A - Dynamic LiNX P326D - PG VR2
Warranty	18 Months : Motor and Electronics 5 year/ limited lifetime frame


OPTIONS

- Unoccupied transit tie down kit
- 22"W x 20" -22"D x 22"H seat
- 7.5" power seat lift
- Swing-away joystick
- USB Charger
- Elevating legrests
- Oxygen tank holder
- Foldable rear basket
- Walker holder
- Cane & crutch holder

MODELS

Vision Sport - Captain	P326-2S C
Vision Sport - Pan	P326-2S S
Vision Sport w/ PSL - Captain	P326D
Vision Sport w/ PSL - Pan	P326D

4 COLOR OPTIONS



Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

Phone: 1-800-963-7487 Fax: 239-772-3252

sales@meritsusa.com

merits