

M A X I C O M F O R T ® S E R I E S



**Introducing the
New MaxiComfort
with a style all
its own!**



Bourbon




Smoke



Brisa Coffee Bean*



Brisa Shiitake*



Brisa Iron*

* upcharge applies

PR735 MaxiComfort
shown in Shiitake

 **Pocket Coil
Spring Seat
Standard!**

EZ^{ZZ}Sleeper

**MAXI
COMFORT®**
Ultimate Recline Technology

**EZ Sleeper
PR735-MLA**

The PR735 EZ Sleeper is just for you! This luxury power lift recliner combines a contemporary design with our patented MaxiComfort positioning for unmatched style and comfort. Its ultramodern, clean lines pair perfectly with a contemporary décor, especially in smaller spaces. Enjoy the long-lasting comfort of this rounded, track arm style chair with a tailored seam back. Your choice of luxurious fabrics and MaxiComfort positioning will make this chair the favorite seat in the house.

 **GOLDEN**

PROUDLY MADE

BY AMERICANS
with Domestic & Imported Materials

Model Number	PR735-MLA
Size	Medium/ Large
Weight Capacity	375 LB
Positioning	MaxiComfort
Twilight Technology	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Bolt
Overall Width	34.5"
Overall Height	45.5"
Overall Depth (Seated)	40"
Width Between Arms	21.5"
Floor to Top of Seat	20"
Seat Depth	22.5"
Seat to Top of Back	30"
Seat to Top of Arm	7.5"
Seat Angle Raised	15°
Distance Required from Wall Reclined	20"
Extension Reclined from Back of Seat to End of Foot Rest	47"
Life Time Warranty, Chair Frame, Lift Frame, Recline Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Warranty on Fabric and Foam, and Pocketed Coils	1 YR
Head Rest Cover	No
Arm Covers	No
Straight Lift	No
Left Hand Control Available	Yes
Extra Pocket Available	Standard
Heat & Massage Available	Yes
Articulating Headrest Available	No
Adjustable Lumbar Available	No
AutoFlex Available	No
Foot Rest Extension Available	Yes
USB Charging Port	No

Battery Backup and Smart Tek Diagnostic System
STANDARD on all UltraComfort models


Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.




Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.

Patented MaxiComfort Positions with Rising Seat



TV



Zero Gravity



Trendelenburg

MaxiComfort patented positioning technology allows for a greater range of positions at the touch of a button. These positions may help improve circulation, relieve back pain and enhance comfort.



shown in Shiitake



Smart Tek

Our luxurious lift and recline chairs are equipped with electronics to help assure a secure lift every time you need it.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with
UL 73 and
CSA C22.2 No. 68

