



### Standard Fabrics:



Calypso



Copper



Hazelnut



Shiraz



Sterling

This fabric stands out for its comfort, smooth texture, and rich color options. It features our exclusive SofTouch™ with Stain Defense coating.

\* An upcharge will apply for all crossover fabrics.

## Relaxer PR766



*Comfort Zone*™  
by GOLDEN

The Relaxer is the original power lift chair recliner featuring our patented MaxiComfort® Ultimate Recline Technology standard! The Relaxer model is a true Golden customer favorite, now with a design refresh with modern, well-appointed features like super plush rounded armrests, an adjustable three-pillow waterfall back, an updated pocket style, and an extra wide chaise pad for better lower leg support. Every MaxiComfort® chair comes standard with our exclusive AutoDrive™ 3.0 hand control to make reclining the chair into your favorite position a breeze!

### Standard Features:

- Patented MaxiComfort® Positioning
- Programmable AutoDrive™ Hand Control
- Extra Pocket



### Quick Ship?

Check real-time inventory at [webportalgt.com](http://webportalgt.com). (Access available for Golden retailers only.)

# MAXICOMFORT® Series PR766

Model Number	PR766-MED
Size	Medium
Weight Capacity	375 LB
Positioning	MaxiComfort
Comfort Zone(s)	2
Twilight Technology	NO
Chaise Seat	YES
Pocketed Coil Seat	YES
Back Type	3 Pillow
Removable Back Type	Bolt
Head Rest Cover	NO
Arm Rest Covers	NO
Overall Width	36"
Overall Height	44.5"
Width Between Arms	20.5"
Floor to Top of Seat	20"
Seat Depth	20.5"
Seat to Top of Back	28.5"
Seat to Top of Arm	5"
Distance Required from Wall Reclined	26"
Extension Reclined from Back of Seat to End of Foot Rest	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Inside Arm
Extra Pocket Option	Standard
Heat & Massage Option	YES
Articulating Headrest Option	NO
Adjustable Lumbar Option	NO
AutoFlex Option	YES
Foot Rest Extension Option	YES
Straight Lift Option	YES
USB Charging Port	NO
Battery Backup	YES
Place of Origin	USA

SmartTek Diagnostic System STANDARD on all Golden models.



**Upgrade with  
Deluxe Therapeutic  
Heat & Massage!**



Shown in  
Hazelnut

## Comfort Zone™ by GOLDEN

Comfort Zones 1 & 2 work together to create individualized wellness positioning and comfort. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.



**WE PATENTED COMFORT™**



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374  
Fax: 800-628-5165  
www.goldentech.com



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with  
UL 73 and  
CSA C22.2 No. 68



www.facebook.com/goldentechnologies

SHT-PR766-031021