

Model	Victory® LX Sport <sup>1</sup>
Product Number	4W: S710LXW
Weight Capacity	400 lbs.
Maximum Speed <sup>2</sup>	Up to 8 mph
Ground Clearance <sup>4</sup>	2.2" at motor, 3.15" at mid-frame
Turning Radius <sup>4</sup>	61.25"
Overall Length <sup>4</sup>	47"
Overall Width <sup>4</sup>	25.5"
Seat-to-Ground Height	23" - 25"
Seat-to-Deck Height	16.5" -18.5"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	Front and rear
Range Per Charge <sup>2,3</sup>	Up to 19.7 miles at 200 lbs. Up to 13 miles at 400 lbs.
Total Weight Without Batteries <sup>6</sup>	160 lbs.
Heaviest Piece When Disassembled <sup>4</sup>	77 lbs. (front section)
Standard Seating	Type: Limited reclining high-back Weight: 43 lbs. Material: Black Vinyl Dimensions: 18" width, 18" depth
Drive System	Rear-wheel drive, sealed transaxle with 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements <sup>5,7</sup>	Size: (2) 50 AH Weight: 30.2 lbs. each
Battery Charger	5 amp, off-board
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months

### HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 4 AGM or gel cell type required.
- 5 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 6 Battery weight may vary based on manufacturer.
- 7 Includes standard seat weight.

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

0001958\_091020



**VICTORY<sup>®</sup>LX<sup>1</sup>**  
**SPORT**

**ENJOY A  
SMOOTH  
AND  
SPORTY  
RIDE**

FDA Class II Medical Device<sup>1</sup>

**LIVE YOUR BEST<sup>®</sup>**  
pridemobility.com |    

# SMOOTH, COMFORTABLE RIDE

18" x 18"  
Vinyl reclining  
high-back seat  
with sliders,  
height adjustable  
headrest and  
flip-up, height  
and width-  
adjustable  
armrests

Puddle light  
under tiller

Standard  
front basket

Easy-to-grip  
tiller adjustment  
knob

Full LED  
lighting  
package  
includes  
curb light,  
rear lights  
and  
directional  
signals

Storage  
compartment  
under the tiller

CTS Suspension includes adjustable  
shocks for greater comfort



Tiller console



CTS Suspension

Maintain your active lifestyle with the **Victory® LX Sport**. Equipped with **Comfort-Trac Suspension (CTS)**, the Victory LX Sport easily handles more varied terrain while providing a smooth and comfortable ride. Enjoy convenient storage and a courtesy light under the tiller. A long-lasting LED headlight ensures optimal pathway illumination. Experience great features and a sporty style with the Victory LX Sport.

- CTS Suspension includes adjustable shocks for greater comfort
- Exclusive one-hand feather touch disassembly
- Charger port located on the tiller lets you conveniently charge your scooter
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrist
- Easy-to-read LED voltmeter
- Rearview mirror

## Available Colors



True Blue



Candy Apple Red

## Optional Seating

Limited recline, black vinyl, high back with height adjustable headrest, flip-up height and width adjustable armrests

- 22" x 18"
- 20" x 20"

## Accessories

- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Additional rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Lap belts: 50", 60" and 70"

## Military Patches

